Hospital for Special Care, New Britain, CT

A Case Study on Ventilator Alarm Management

Hospital for Special Care

Presenter: Connie Dills, MBA, RRT, RPFT

Hospital for Special Care

is 228 bed Long Term Acute Care Hospital (LTACH)

We Rebuild Lives.



11 bed regional ventilator weaning unit

Hospital for Special Care

A Case Study on Ventilator Alarm Management

is 228 bed Long Term Acute Care Hospital (LTACH)

We Rebuild Lives.



27 bed pediatric unit

Hospital for Special Care

A Case Study on Ventilator Alarm Management

is 228 bed Long Term Acute Care Hospital (LTACH)

We Rebuild Lives.



2 chronic ventilator units (36 & 38 beds)

Hospital for Special Care

A Case Study on Ventilator Alarm Management

is 228 bed Long Term Acute Care Hospital (LTACH)

We Rebuild Lives.



10 bed cardiac unit -CHF & VADs

Hospital for Special Care

A Case Study on Ventilator Alarm Management

is 228 bed Long Term Acute Care Hospital (LTACH)

We Rebuild Lives.



28 bed rehabilitation unit -CVA, pulmonary, TBI, SCI, Orthopedic

> Hospital *for Special* Care

is 228 bed Long Term Acute Care Hospital (LTACH)

We Rebuild Lives.



23 bed satellite unit in Hartford

Hospital for Special Care

A Case Study on Ventilator Alarm Management

is 228 bed Long Term Acute Care Hospital (LTACH)

We Rebuild Lives.



8 bed autism unit

Hospital *for Special* Care

A Case Study on Ventilator Alarm Management

is 228 bed Long Term Acute Care Hospital (LTACH)

We Rebuild Lives.



2 neurobehavioral units (10 & 15 beds)

Hospital *for Special* Care

A Case Study on Ventilator Alarm Management

is 228 bed Long Term Acute Care Hospital (LTACH)

We Rebuild Lives.



22 bed chronic rehabilitation unit

Hospital for Special Care

A Case Study on Ventilator Alarm Management

The Issue: Managing 100 ventilators a day

- A typical critical care ventilator can alarm for approximately 135 different reasons
- Ventilator population at HSC comprised primarily of awake & active patients, most of which are in private & semi-private rooms
- Average number of ventilator device alarms per ventilator per day at HSC is a staggering 48 X 100 = 4800!

Hospital for Special Care

The Solution

Middleware

2003 the organization sought options for increasing patient safety by utilizing an alarm management system

- Existing monitoring was through nurse call
- Every ventilator alarm, every time
- Bernoulli® was the product selected
- Alarms filtered (LowVe, LIP, No Data, Patient Disconnect)
- Immediately actionable
- Not immediately actionable but if unresolved default to actionable



perturbation of the second

How We Manage: The Basics

- The Respiratory Therapists use Bernoulli® for their patient assignment
- Each patient is assigned to the Respiratory Therapist taking care of them for that shift by assigning a pager that the RT carries
- When actionable alarm conditions occur the RT is alerted by pager, overhead audible alarm, visual alerts through desktop/laptop computers and in pediatrics via an overhead scrolling message bar

Alarm settings are verified by RTs once a shift

Hospital for Special Care

A Case Study on Ventilator Alarm Management











Enhanced ventilator alarm response time <u>18 second average</u>

Hospital for Special Care

Slide 20

A Case Study on Ventilator Alarm Management

Realized an <u>80% reduction</u> in alarms requiring immediate action

Hospital for Special Care

Slide 21

A Case Study on Ventilator Alarm Management

No serious ventilator related safety events in 13 years

Hospital for Special Care

Slide 22

A Case Study on Ventilator Alarm Management

Significant reduction in ambient noise and resultant alarm fatigue

Hospital for Special Care

Slide 23

A Case Study on Ventilator Alarm Management

The Joint Commission NPSG.06.01.01 Improve the Safety of Clinical Alarm Systems



A Case Study on Ventilator Alarm Management

NPSG.06.01.01: Alarm Management

- Monitoring and responding to alarm signals
- Checking individual alarm signals for accurate settings, proper operation, and detectability
- Customize for specific clinical units, groups of patients, or individual patients

Joint Commission Perspectives®, July 2013, Volume 33, Issue 7 Copyright 2013 The Joint Commission

> Hospital for Special Care

Current Conditions: Well Positioned

In meeting TJC's NPSG.06.01.01 by 1/1/2016

Middleware facilitates monitoring & checking of alarm signals for accurate settings, proper operation & detectability for appropriate response on a shift to shift basis

> Hospital for Special Care

A Case Study on Ventilator Alarm Management

The Future: Enhanced Patient Safety

Data Gathering & Benchmarking

Response time

Reducing alarm duration - fatigue

Further alarm reduction efforts

- Device specific bedside for patient comfort
 - Enhanced ventilator technology
 - Device alarm volumes

Benchmarking with other sites

Hospital for Special Care