

Nurses and Respiratory Therapists – Working Together for Safe Alarm Systems Management

May 11, 2015

Speaker Introductions

- **Shawna Strickland**, PhD, RT-NPS, RRT-ACS, AE-C, FAARC, Associate Executive Director Education, American Association for Respiratory Care
- **Jenifer L. Burke**, RN, MSN, CPNP, Pediatric Pulmonary & Cystic Fibrosis, Rush University Medical Center

Disclosures

Jenifer Burke

- Employee, Rush University Medical Center
- Adjunct Faculty, College of Nursing, Graduate College, Rush University

Shawna Strickland

- Employee, American Association for Respiratory Care
- Contract, Centers for Disease Control and Prevention, SNS workshop
- Adjunct Faculty, College of Health Sciences and Graduate College, Rush University

Objectives

- Describe the basics of ventilator alarms
- Rank the importance of ventilator alarms based on the physiologic responses the ventilator alarms represent
- Develop a strategy to respond appropriately to various ventilator alarms
- Identify the daily nursing activities that trigger ventilator alarms

Alarm Safety and Fatigue

The Problem

- Visual/audible alarms when the patient's condition changes or machine error
- Drastically increased number of devices with audible alarms at the bedside.
- Overwhelmed bedside practitioners exposed to different levels of audible alarms

The Effects

- Alarm fatigue
- “False alarms”
- Nuisance alarms
- Adverse patient outcomes
- #1 on ECRI Institute Top 10 Health Technology Hazards in 2015
- TJC National Patient Safety Goals

Mechanical Ventilation 101

Indications

- Apnea
- Acute respiratory failure
- Impending respiratory failure
- Refractory hypoxemia

Types

- Invasive
 - Requires artificial airway (endotracheal tube or tracheostomy tube)
 - Provides airway for patient who cannot protect his/her own airway
- Non-invasive
 - Delivered via face mask, nasal mask/pillows
 - Does not provide a protected airway
 - Patient must be able to protect his/her own airway

Mechanical Ventilation 101

Pressure v Volume

- Pressure:
 - Breath terminates when pre-set pressure is reached
 - Volume is variable depending on patient compliance and resistance
- Volume
 - Breath terminates when pre-set volume is reached
 - Pressure is variable depending on patient compliance and resistance

Modes

- Pressure
 - PC-CMV, PC-SIMV, PSV, AVAPS
- Volume
 - VC-CMV, VC-SMIV
- Dual modes
 - Pressure limited, volume targeted (VS, PRVC)
 - Pressure limited, volume guaranteed


Mechanical Ventilation 101

Breath Types

- Spontaneous
 - Patient initiates, patient determines depth and length
- Supported
 - Patient initiates, machine supports depth
- Mandatory
 - Machine initiates, machine determines depth and length


Trigger Types

- Pressure
 - Preset pressure detected
- Flow
 - Preset flow detected
- Volume
 - Preset volume detected
- Time
 - Preset time interval has elapsed




Patient effort;
patient triggers
machine with
negative
pressure


No patient
effort; no
negative
deflection below
pressure baseline



Pressure Support Ventilation



- A. Patient effort
- B. Support from machine (PS)
 - B1-Over-shoot
 - B2-Under-shoot
- C. Plateau
- D. Termination of support



Mechanical Ventilation 101

Pressure Settings

- Respiratory rate (f)
- Peak pressure (PIP)
- Inspiratory time (T_I)
- Positive expiratory pressure (PEEP)
- Fraction of inspired oxygen (FiO_2)



Volume Settings

- Respiratory Rate (f)
- Tidal volume (V_T)
- Inspiratory flow (V)
- Positive expiratory pressure (PEEP)
- Fraction of inspired oxygen (FiO_2)

Mechanical Ventilation 101


- Measured values
 - Peak inspiratory pressure (PIP)
 - Plateau pressure (P_{PLAT})
 - Minute ventilation (V_E)
 - Auto-PEEP
 - Total respiratory rate
 - Exhaled tidal volume (V_T)

Anatomy of a Waveform




Influencing Factors


Oxygen



Carbon Dioxide



Patient in Distress



Potential Ventilator Alarms

- High pressure
 - Achieved PIP is too high
- Low pressure
 - Achieved PIP is too low
- High PEEP
 - Measured PIP is too high
- Low PEEP
 - Measured PIP is too low
- Apnea
 - RR falls below set threshold
- Inverse I:E ratio
 - Inspiration is longer than exhalation
- High tidal volume
 - Exhaled V_T is too high
- Low tidal volume
 - Exhaled V_T is too low
- High minute volume
 - Exhaled V_E is too high
- Low minute volume
 - Exhaled V_E is too low
- High/low respiratory rate
 - Patient total RR too high/low
 - Includes spontaneous rates


What did I do?

↑ V_T , V_E or RR alarm

- Air hungry
- Sigh
- Pain
- Agitation
- Under sedation
- Procedures
- Water in tube

What do I do?

↑ V_T , V_E or RR alarm




What did I do?

↓ pressure, PEEP, VT, V_E alarm

- Disconnected the vent
- Didn't inflate cuff
- Suctioning
- Over sedation (spontaneous modes)
- Leak in circuit


What do I do?

↓ pressure, PEEP, VT, V_E alarm




↓ pressure, PEEP, VT, V_E alarm: Air Leaks

Volume-Time Scalar



Pressure-Volume Loop




What did I do?

↑ pressure or PEEP alarm

- Patency of tube
(blocked/clamped/bent)
- Secretions
- Cough
- Resistance
- Poor positioning


What do I do?

↑ pressure or PEEP alarm



What do I do?


↑ pressure or PEEP alarm




↑ pressure or PEEP alarm

worsening C_L

Pressure Ventilation



Volume Ventilation




Copyright © 2004, 1999, Mosby, Inc. All Rights Reserved.


Copyright © 2004, 1999, Mosby, Inc. All Rights Reserved.

↑ pressure or PEEP alarm worsening R_{AW}

Much longer
expiratory time;
increased R_{AW} =
takes longer to
exhale




A



Area within
the loop
(hysteresis)
is much
larger

C



B


Reduced
expiratory
flow
(speed of
exhalation)

(From Pilbeam SP: *Mechanical ventilation: physiological and clinical applications*, ed 3, St Louis, 1998, Mosby.)

Copyright © 2004, 1999, Mosby, Inc. All Rights Reserved.

What do I do?

↑ pressure or PEEP alarm




What did I do?

Inverse I:E Ratio alarm

- Anything that would cause a change in the RR
- Patient has an increased drive to breathe
 - Pain
 - Need for sedation

What do I do?

Inverse I:E Ratio alarm




What did I do?

Apnea alarm

- Over sedated
- Disconnected
- Patency of tube


What do I do? Apnea alarm



Apnea alarm


Sensitivity

Pressure drops below baseline because patient has to pull to trigger breath



Air Leaks

Exhaled volume does not reach baseline



AAM

Case Study

- 23 year old male
- Acute asthma exacerbation secondary to viral infection
- Intubated in ED with 7.5 ETT
- Currently receiving:
 - Midazolam
 - Fentanyl
 - Albuterol via nebulizer
- Ventilator settings:
 - VT = 600 mL (6ml/kg PBW)
 - RR = 10 breaths/min
 - FiO₂ = 1.0 (100%)
 - PEEP = 0
 - Flow = 90 L/min



Case Study: After Initiation of Mechanical Ventilation

- Blood gases:

- pH = 7.35
- PaCO₂ = 47 mm Hg
- PaO₂ = 186 mm Hg
- HCO₃ = 25 mEq/L
- SaO₂ = 91%

- Measured Ventilator Parameters:

- PIP = 65 cm H₂O
- Plateau = 25 cm H₂O
- AutoPEEP = 8 cm H₂O



Copyright © 2006, Mosby Inc.

(From Dhand R: Respir Care 50:256, 2005.)

Take-Home Points

- Always look at the patient first
- Don't fall victim to alarm fatigue
 - Know your ventilators (capabilities, sounds, etc.)
 - Know your established protocols
 - Confirm ventilator settings and alarm settings at the when you first accept the patient for your shift
- When in doubt and the patient is in distress, manually ventilate the patient and call for help

References

- Bridi AC, Louri TQ, da Silva RC. Clinical alarms in intensive care: implications of alarm fatigue for the safety of patients. Rev Lat Am Enfermagem 2014;22(6):1034-1040.
- Cairo JM, Pilbeam SP. Mosby's respiratory care equipment, 7th edition. St. Louis: Mosby|Elsevier; 2004.
- Dunlevy C, Strickland S. Respiratory therapists' response to a national clinical alarm survey. AARC Times 2012, 42-48.
- ECRI Institute. Top 10 health technology hazards for 2015. Health Devices 2014.
- Funk M, Clark T, Bauld TJ, Ott JC, Coss P. Attitudes and practice related to clinical alarms. Am J Crit Care 2014;23(3):e9-e18.
- Hess DR, McIntyre NR, Mishoe SC, Galvin WF, Adams AB, Saposnick AB. Respiratory care principles & practices. Elsevier; 2002.

References

- Hess DR, McIntyre NR, Mishoe SC, Galvin WF, Adams AB, Saposnick AB. Respiratory care principles & practices. Elsevier; 2002.
- Joint Commission on Accreditation of Healthcare Organizations. Sound the alarm: managing physiologic monitoring systems. The Joint Commission Perspectives on Patient Safety 2011;11(12):6-11.
- Lowery WS. Ventilator-disconnect and death: a case study and a safety device. Respir Care 2010;55(6):774-776.
- Papadakos PJ. Electronic distractions of the respiratory therapist and their impact on patient safety. Respir Care 2014;59(8):1306-1309.
- Pilbeam SP. Mechanical ventilation: physiological and clinical applications, 4th edition. Mosby|Elsevier; 2006.
- Walsh BW, Czervinske MP, DiBlasi RM. Perinatal and pediatric respiratory care, 3rd edition. Elsevier; 2010.

Thank you!

Jenifer Burke RN, MSN, CPNP

Pediatric Nurse Practitioner
Rush University Medical Center

Jennifer_Burke@rush.edu

Shawna Strickland PhD RRT-NPS RRT-ACCS AE-C FAARC

Associate Executive Director-Education
American Association for Respiratory Care
Adjunct Faculty, College of Health Sciences and Graduate
College
Rush University

shawna.strickland@aarc.org

AAMI FOUNDATION HTSI

Healthcare Technology Safety Institute

Thank You to Our Industry Partners

This Patient Safety Seminar is offered at no charge thanks to commercial grants from the companies shown here. The AAMI Foundation and its co-convening organizations appreciate their generosity. The AAMI Foundation is managing all costs for the series. The companies had no role in content development, and the seminars do not contain commercial content. Seminar presenters were selected based on topic expertise without regard to industry affiliation by a multi-disciplinary AAMI Foundation volunteer planning committee.

PLATINUM



COVIDIEN



GE Healthcare

PHILIPS

GOLD

Baxter



Dräger

mindray

SILVER



EarlySense

